


UNIQUE FEATURES OF THE *invertrac*®

- Easy to mount and dismount. The thigh pad is attached to a crossbar which has a low center of gravity allowing the user complete control when inverting or dismounting. This unique design makes this an ideal unit for patient in-home use.
- The thigh pad has been designed to reduce pressure in the pubic area during inversion sessions and when performing extension/flexion exercises.
- Rubber sidegrips for safe and easy inversion and dismount.
- Fully adjustable for patient leg length, thigh thickness and height - maximizing comfort and traction during each inversion session.
- M-Beam at front allows for extra-gentle extensions.
- Designed and built for the professional. Chrome-plated heavy duty steel tubing and vinyl components will give years of trouble free use.
- 99% assembled - ready to use just minutes from unpacking. Simply attach the thigh pad and slip on the adjustable restraining bar.
- Folds completely flat for easy transportation and storage.
- Lightweight, approximately 33 pounds.
- Tested to safely hold up to 350 pounds.
- Accommodates most individuals up to 6'4". Optional kit extends to 6'8".

STEP 1




EASY TO INVERT...


◀ Step into unit and ensure pelvic (hip) bones are slightly above the top of thigh support/seat. Height adjustment may be required.

▶ Lean forward while holding handlegrips. To raise restraint bar, raise one foot.

STEP 2




STEP 3



◀ Grasp side-grips and tip forward. Continue to M-Beam until fully inverted. Relax...

▶ This is the correct position for natural traction. Head should not be touching floor, place hands in a relaxed position, as shown. Adjust height and femur dial if necessary.

STEP 4



TESTIMONIALS

"Convenient device for performing a full range extension exercise... leads to improvement in back function."

Murray E. Allen M.D.

Associate Prof. Kinesiology
Simon Fraser University

★ ★ ★

"The Invertrac offers many of our clients relief from problematic compression injuries, and often provides a more effective forward flexion stretch than stretching alone."

Kim N. Zimmer B.Sc. Kin

Kinesiologist
Metro Rehab Centre

★ ★ ★

"I have found the Invertrac to be the superior product. The Invertrac is definitely less expensive, easier for the patients to mount and dismount, easier to store and carry, and more comfortable to use."

R.H. Latch, B.Sc., D.C.

Inglewood Chiropractic Office

★ ★ ★

"Lightweight and quite portable, I have used it with athletes and coaches alike... very beneficial where traction is indicated."

Carl Petersen, Physiotherapist

National Alpine Ski Team

★ ★ ★

"Very effective in reducing chronic low back conditions... frequency of recurring problems is often reduced."

Dr. I. Gainor,

Calgary, AB

★ ★ ★

"Due to your unit's easy inversion and dismount features I have been able to utilize the Invertrac with patients who have undergone spinal surgery, resulting in a much faster recovery period."

Dr. M. Shales,

Alabama

Distributor:

SYNERGEXX HOLDINGS INC.

Canada

Toll Free:

1-800-667-7795

Website:

www.invertrac.com

Email:

invertrac@invertrac.com

Tel/Fax:

On website